

pSG5 Vector

INSTRUCTION MANUAL

Catalog #216201

Revision #111001a

For In Vitro Use Only



* 2 1 6 2 0 1 - 1 2 1 1 1 0 0 1 A / *

LIMITED PRODUCT WARRANTY

This warranty limits our liability to replacement of this product. No other warranties of any kind, express or implied, including without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided by Stratagene. Stratagene shall have no liability for any direct, indirect, consequential, or incidental damages arising out of the use, the results of use, or the inability to use this product.

ORDERING INFORMATION AND TECHNICAL SERVICES

United States and Canada

Stratagene

11011 North Torrey Pines Road

La Jolla, CA 92037

Telephone (858) 535-5400**Order Toll Free** (800) 424-5444**Technical Services** (800) 894-1304**Internet** tech_services@stratagene.com**World Wide Web** www.stratagene.com

Stratagene European Contacts

Location	Telephone	Fax	Technical Services
Austria	0800 312 526	0800 312 527	017 956 7036
Belgium	00800 7000 7000	00800 7001 7001	00800 7400 7400
	0800 96078	0800 96024	027 13 12 11
France	00800 7000 7000	00800 7001 7001	00800 7400 7400
	0800 100391	0800 881323	
Germany	00800 7000 7000	00800 7001 7001	00800 7400 7400
	06221 400 634	06221 400 639	06221 400 637
	0130 840 911	0130 762 088	069 9509 6197
Netherlands	00800 7000 7000	00800 7001 7001	00800 7400 7400
	0800 023 0446	0800 023 0447	0800 023 0448
Switzerland	00800 7000 7000	00800 7001 7001	00800 7400 7400
	061 6930 540	061 6930 544	01 800 9045
	0800 830 250	0800 825 225	
United Kingdom	00800 7000 7000	00800 7001 7001	00800 7400 7400
	01223 420 955	01223 420 234	0171 365 1056
	0800 585 370	0800 783 0889	

All Other Countries

Please contact your local distributor. A complete list of distributors is available at www.stratagene.com.

pSG5 Vector

CONTENTS

Materials Provided	1
Storage Conditions	1
Vector Features	2
Applications.....	2
The pSG5 Vector	3
Preparation of Host Strain	4
Preparation of a –80°C Glycerol Stock	4
Preparation of Media and Reagents	4

pSG5 Vector

MATERIALS PROVIDED

Material Provided	Quantity
pSG5 vector	20 μ g
AG1 strain*, glycerol stock	1 tube

* *recA1*, *endA1*, *gyrA96*, *thi-1*, *hsdR17*, (*r_k⁻*, *m_k⁺*), *supE44*, *relA1*, (uncharacterized mutation improves transformation efficiency)

STORAGE CONDITIONS

pSG5 vector: -20°C

AG1 bacterial glycerol stock: -80°C

Revision #111001a

Copyright © 2001 by Stratagene

VECTOR FEATURES


The pSG5 Vector is a eukaryotic expression vector constructed by combining pKCR2 and Stratagene's pBS vector. Because of the high copy number of this plasmid, large quantities of double-stranded DNA are obtained.

Applications

The pSG5 is useful for both in vitro and in vivo expression. Expression in vivo is achieved via transient expression in a variety of cell lines (highest level of expression is obtained following transfection of a cell line expressing the T antigen). The f1 origin allows rescue of ssDNA for use in mutagenesis and sequencing. The SV40 early promoter and polyadenylation signal promotes expression in vivo, and the T7 bacteriophage promoter facilitates in vitro transcription of cloned inserts. The β -globin intron II allows splicing of expressed transcripts.

To ligate the gene of interest into the pSG5 vector, use the unique restriction sites *EcoR* I, *Bam*H I, and *Bgl* II (downstream from the promoter).

The pSG5 Vector



Feature	Nucleotide Position
SV40 promoter and SV40 origin of replication	28–366
β-globin intron	395–967
T7 promoter	1022–1040
EcoR I	1043
BamH I	1049
Bgl II	1055
SV40 polyA signal	1069–1202
pUC origin of replication	1342–2009
ampicillin resistance (<i>bla</i>) ORF	2160–3017
f1 origin of ss-DNA replication	3587–3893

Figure 1 Circular map and features of the pSG5 vector. The complete sequence and list of restriction sites is available at www.stratagene.com.

PREPARATION OF HOST STRAIN

The host strain has been sent as a glycerol stock. For the appropriate media and plates, please refer to the following table:

Bacterial strain	Plates for bacterial streak	Media for glycerol stock
AG-1	LB	LB

On arrival, prepare the following from the glycerol stock:

Note *Do not allow the contents of the vial to thaw. The vials can be stored at -20 or -80°C , but most strains remain viable longer if stored at -80°C .*

1. Revive the stored cells by scraping off splinters of solid ice with a sterile wire loop.
2. Streak the splinters onto an LB plate.

Restreak the cells fresh each week.

Preparation of a -80°C Glycerol Stock

1. In a sterile 50-ml conical tube, inoculate 10 ml of the appropriate liquid media with one or two colonies from the plate. Grow the cells to late log phase.
2. Add 4.5 ml of a sterile glycerol–liquid media solution (5 ml of glycerol + 5 ml of liquid media) to the bacterial culture from step 1. Mix well.
3. Aliquot into sterile centrifuge tubes (1 ml/ tube).

This preparation may be stored at -20°C for 1–2 years or at -80°C for more than 2 years.

PREPARATION OF MEDIA AND REAGENTS

LB Broth (per Liter)

10 g of NaCl
10 g of tryptone
5 g of yeast extract
Add deionized H_2O to a final volume of
1 liter
Adjust to pH 7.0 with 5 N NaOH
Autoclave

LB Agar (per Liter)

10 g of NaCl
10 g of tryptone
5 g of yeast extract
20 g of agar
Add deionized H_2O to a final volume of 1 liter
Adjust pH to 7.0 with 5 N NaOH
Autoclave
Pour into petri dishes (~25 ml/100-mm plate)