



NYCOMPS expression vector

1.

CATATCCTGCAGGCGCGCCAAATTGCGGCCGCGGCCGAGGAAAACCTGTATTTTCAGGGATCCCATCATCAACCATCACCATCACCATCACTAACTCGAG

A A A A G Q E N L Y F Q G S H H H H H H H H H *

NdeI SbfI AscI NotI FseI TEV protease BamHI His*10 STOP XhoI

GENE of interest EGFP

This linker will be inserted into pET24, which is kanamycin resistant, using the NdeI and XhoI sites. EGFP will insert, in frame at the BamHI site. The gene of interest will be inserted into the vector using SbfI and FseI or AscI and NotI double digests. The ATG of the gene of interest will determine the frame. The stop codon should be removed from the gene of interest. This vector may also be used, without the EGFP inserted, as a His- tagged vector only.